# qUAntify<sup>®</sup> Plus Control Levels 1 and 2

REF 995 Bilevel 10 x 12 mL  
995X Bilevel MiniPak 2 x 12 mL

EXP 2022-06-03

LOT 80670 Level 1 80671  
Level 2 80672

## ENGLISH

### INTENDED USE

qUAntify<sup>®</sup> Plus Control is intended for use as an assayed quality control to monitor the precision of urinalysis test procedures for the analytes listed in this package insert.

### SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. Two levels of control are available to allow monitoring of performance.

For customers in Germany: Quality control materials are required for assessment of laboratory performance as described in the "Guideline for Quality Assurance of Medical Laboratory Examinations following the German Medical Association" (Rili-BÄK regulation).

### REAGENT

This product is a liquid base matrix prepared with human erythrocytes and leukocytes, constituents of animal origin, chemicals, and preservatives. The control is provided in liquid form for convenience.

### STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened at 2 to 8 °C.

#### Chemical Analysis with Multiple Immersion Use

Once opened, this product will be stable for 14 days at 2 to 8 °C or 8 test strip immersions, whichever occurs first.

#### Microscopic and Chemical Analysis / Pregnancy (hCG) Testing (Non-centrifuged)

Once opened, this product will be stable for 31 days when stored tightly capped at 2 to 25 °C with the following exception: Ketones will be stable for 31 days at 2 to 8 °C or 10 days at room temperature (18 to 25 °C). For optimum performance, store the product at 2 to 8 °C.

### DO NOT FREEZE. Store this product away from the light.

Note: It is normal to observe sediment constituents in the bottom of the tube when stored for a prolonged period.

This product is shipped under refrigerated conditions.

### PROCEDURE

This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit, or reagent being used.

Before sampling, allow the tube to reach room temperature (18 to 25 °C) for 15 to 30 minutes.

Invert the tube several times to ensure homogeneity. **DO NOT SHAKE.**

#### Chemical Analysis with Multiple Immersion Use

1. Remove cap and immerse the urinalysis test strip in the tube.
2. Read the test strip, visually or with an instrument reader, in accordance with manufacturer's instructions.
3. After each use, promptly recap the tube and return it to refrigerated storage conditions.
4. Repeat steps 1 to 3 for a total number of no more than 8 immersions.

#### Microscopic and Chemical Analysis (one time use)

1. Remove cap and immerse the urinalysis test strip in the tube.
2. Read the test strip, visually or with an instrument reader, in accordance with manufacturer's instructions.
3. Perform microscopic analysis in accordance with established laboratory procedures.
4. Discard the remaining material.

#### Pregnancy (hCG) Testing

After each use, promptly recap the tube and return it to refrigerated storage conditions.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.

### LIMITATIONS

1. This product should not be used past the expiration date.
2. If there is evidence of microbial contamination or excessive turbidity in the product, discard the controls.
3. This product is not intended for use as a standard.
4. Values may gradually decrease over the product shelf life for the following analyte: Ketones
5. Values may gradually increase over the product shelf life for the following analyte: pH

### ASSIGNMENT OF VALUES

The results printed in this insert were derived from replicate analyses and are specific for this lot of product. The tests listed were performed by the manufacturer and/or independent laboratories using manufacturer supported reagents and a representative sampling of this lot of control. Variations over time and between laboratories may be caused by differences in laboratory technique, instrumentation, and reagents, or by manufacturer test method modifications. Each laboratory should use the results provided only as a reference and establish its own parameters of precision.

Refer to [www.qcnet.com](http://www.qcnet.com) for insert update information.

### SPECIFIC PERFORMANCE CHARACTERISTICS

This product is a stabilized liquid product manufactured under rigid quality control standards. To obtain consistent assay values, the control requires proper storage and handling as described.

### OPTIONAL SUPPLIES

Data Log Sheets available through local subsidiary or at [www.qcnet.com](http://www.qcnet.com).

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Analyte	Manufacturer/Method	Units	Level 1 – 80671	Level 2 – 80672	Units	Level 1 – 80671	Level 2 – 80672
<b>RED BLOOD CELLS (RBC)</b>							
ARKRAY AUTION Hybrid AU-4050		Cells/hpf (5)	0 - 2	5 - 40	Cells/ $\mu$ L	0 - 10	50 - 200
Iris Diagnostics IQ200 Series		Cells/hpf (5)	0 - 5	5 - 35	Cells/ $\mu$ L	0 - 28	28 - 193
Sysmex UF-50/100/1000		Cells/hpf (5)	0 - 5	5 - 35	Cells/ $\mu$ L	0 - 28	28 - 194
KOVA System (6)		Cells/hpf (5)	0 - 5	0 - 60	Cells/hpf (5)	0 - 5	0 - 60
Sediment / Slide & Coverslip		Cells/hpf (5)	0 - 5	0 - 60	Cells/hpf (5)	0 - 5	0 - 60
<b>WHITE BLOOD CELLS (WBC)</b>							
ARKRAY AUTION Hybrid AU-4050		Cells/hpf (5)	0 - 2	5 - 30	Cells/ $\mu$ L	0 - 10	50 - 200
Iris Diagnostics IQ200 Series		Cells/hpf (5)	0 - 5	5 - 60	Cells/ $\mu$ L	0 - 28	28 - 330
Sysmex UF-50/100/1000		Cells/hpf (5)	0 - 5	0 - 30	Cells/ $\mu$ L	0 - 28	0 - 167
KOVA System (6)		Cells/hpf (5)	0 - 5	0 - 60	Cells/hpf (5)	0 - 5	0 - 60
Sediment / Slide & Coverslip		Cells/hpf (5)	0 - 5	0 - 55	Cells/hpf (5)	0 - 5	0 - 55
<b>CRYSTALS (7)</b>							
ARKRAY AUTION Hybrid AU-4050		lpf (5)	Absent	Absent	lpf (5)	Absent	Absent
Iris Diagnostics IQ200 Series		hpf (5)	Absent	Absent - Occasionally Present	hpf (5)	Absent	Absent - Occasionally Present
Sysmex UF-50/100/1000		hpf (5)	Absent - Occasionally Present	Absent - Occasionally Present	hpf (5)	Absent - Occasionally Present	Absent - Occasionally Present
KOVA System (6)		hpf (5)	Absent	Absent - Occasionally Present	hpf (5)	Absent	Absent - Occasionally Present
Sediment / Slide & Coverslip		hpf (5)	Absent	Absent - Occasionally Present	hpf (5)	Absent	Absent - Occasionally Present
<b>CASTS</b>							
ARKRAY AUTION Hybrid AU-4050		lpf (5)	Absent	Absent	lpf (5)	Absent	Absent
Iris Diagnostics IQ200 Series		lpf (5)	Absent	Absent	lpf (5)	Absent	Absent
Sysmex UF-50/100/1000		lpf (5)	Absent - Occasionally Present	Absent - Occasionally Present	lpf (5)	Absent - Occasionally Present	Absent - Occasionally Present
KOVA System (6)		lpf (5)	Absent	Absent	lpf (5)	Absent	Absent
Sediment / Slide & Coverslip		lpf (5)	Absent	Absent	lpf (5)	Absent	Absent

## PREGNANCY

Schwangerschaft // Grossesse // Gravidanza // Embarazo // Gravidez // Graviditet // Graviditet // Gebelik // 妊娠

Analyte	Manufacturer/Method	Level 1 – 80671	Level 2 – 80672
<b>hCG</b>			
Alere (Clearview) hCG COMBO II		Negative	Positive
Polymedco Poly stat hCG (HCG25)		Negative	Positive
Quidel QuickVue One Step hCG Combo Test		Negative	Positive
Siemens Clinitest hCG (Siemens Clinitek Status / Status+)		Negative	Positive
Other (8)		Negative	Positive

## FOOTNOTES // Fussnoten // Notes de bas de page // Note a piè pagina // Notas al pie de página // Notas de rodapé // Fotnoter // Fodnoter // Dipnotlar // 脚注

### ENGLISH

- Tested and found to be negative.
- Colors produced by the Urobilinogen and/or Bilirubin test pads may not be characteristic of those shown on the manufacturer's color chart or label. Intensity of color may be correlated to concentration.
- Values apply to all Chemstrip Reagent Strips.
- Values apply to all Siemens Reagent Strips.
- lpf = low power field. hpf = high power field.
- Results were obtained by following the manufacturer's instructions for the KOVA System for Standardized Urinalysis.
- Crystals may appear atypical in size and shape.
- Other test kits with sensitivities of ≥ 10 mIU/mL.
- The data required to establish the means and acceptable ranges for this assay were not obtained due to limited assignment participation. If your facility is interested in participating in the Value Assignment Program for this assay, please contact your local Bio-Rad office.
- INTERNATIONAL USE ONLY** - The following section contains data for methods that are not available for diagnostic use in the United States.

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