Sorval ST-16 Centrifuge Timer Verification Job Aid

| Step | Action | | | | | | |
|--|--|--|--|--|--|--|--|
| **Document all maintenance performed on Centrifuge Maintenance Log** | | | | | | | |
| 1 | Obtain a NIST calibrated/certified stopwatch/timer | | | | | | |
| 2 | Prepare the centrifuge for normal operation, close the lid and select START | | | | | | |
| 3 | Start the stopwatch/timer | | | | | | |
| 4 | Stop the stopwatch/timer immediately when the centrifuge begins deceleration | | | | | | |
| 5 | Compare the set centrifugation time with the stopwatch/timer; document the set time and stopwatch/timer reading on the Centrifuge Maintenance Log and indicate Pass or Fail. Stopwatch/timer must be within ±5% of the set time to pass | | | | | | |
| 6 | Repeat steps 1-5 for all programs if times are different, documenting all set times and measure times as well as Pass or Fail for each measure time | | | | | | |
| 7 | If centrifuge fails timer check, call service and document actions taken on Centrifuge Maintenance Log | | | | | | |

Example of Centrifuge Maintenance Log Sheet

| Semi-Annual – Rpm & Tir RPM Setting | | rpm ± 10% | | ss/Fail | rpm ± 10% | Pass/Fail | Performed By (Initial): | Date |
|--|-----|------------|---|---------|------------|-----------|-------------------------|------|
| Tachometer Reading | | | | | | | | |
| Tachometer Reading | | | | | | | | |
| Timer Check Setting | Mir | Minutes±5% | | s/Fail | Minutes±5% | Pass/Fail | Performed by (initial): | Date |
| Stopwatch reading | 4 | , | 4 | • | | | | |
| Stopwatch reading | | | | | | | | |
| | | | | | | | • | |
| | | | | | | | | |
| | | | | | | | | |

Document set time and actual readings and if Pass or Fail (use additional columns if needed)



JA270-30-01E V01 Effective: 27-NOV-2020 Approved by: H Malvern