

Document History:

Title: Validation of New QC and Strip Lots for Nova StatStrip Glucose Meter

Site(s):

Quality and Development Technologist

| | | | |
|--------------------|----------------------------------|--------------------|---------------|
| Document #: | 110-130-09 | Version #: | 03 |
| Section: | Clinical Biochemistry & Genetics | Subsection: | Point of Care |

| | | | |
|---------------------|--------------|--------------|-------------|
| Approved by: | L Thorlacius | Date: | 24-May-2018 |
|---------------------|--------------|--------------|-------------|

| | | | |
|---------------------|----------------------|------------------------|-------------|
| Approved by: | L Thorlacius | Date: | 29 Jan 2021 |
| Signature: | (Signature on File)_ | Effective Date: | 29 Jan 2021 |

| # | Details of Revisions: | Approval: | Date: |
|---|--|--------------|-------------|
| 1 | New document | L Thorlacius | 24 May 2018 |
| 2 | Updated References to include SOP 110-130-07. Removed Roche Result 2 and Mean from Patient Comparison chart in Appendix C (p.6). | A Sokoro | 26 Jan 2021 |
| 3 | | | |
| 4 | | | |

DISCLAIMER: Please be advised that printed versions of any policy, or policies posted on external web pages, may not be the most current version of the policy. Although we make every effort to ensure that all information is accurate and complete, policies are regularly under review and in the process of being amended and we cannot guarantee the accuracy of printed policies or policies on external web pages. At any given time the most current version of any Shared Health Inc. policy will be deemed to apply. Users should verify that any policy is the most current policy before acting on it.

1.0 REFERENCES:

- 1.1 StatStrip Blood Glucose POC Meter Procedure Manual, SOP # 110-130-01
- 1.2 Evaluation Protocol of POC Blood Glucose Monitoring Instruments, SOP # 110-130-07
- 1.3 Nova Biomedical StatStrip Express Glu/Ket Control 1/3 package insert, Ref# 46947/46949 2017-03
- 1.4 EP Evaluator Software 11.1.0.26, 4161-A6-7B Guide

2.0 PURPOSE:

New lot numbers of Nova quality control material and reagent test strips must be validated in-house prior to provincial distribution.

Nova Biomedical will contact us well in advance indicating there are either new lots of QC or test strips that require validation. The new products are to be sent directly to the Quality and Development Technologist for evaluation. Once validated, Nova Biomedical is notified the product is acceptable for patient testing.

3.0 MATERIALS:

| Reagents | Supplies | Equipment |
|--|--|--|
| <ul style="list-style-type: none"> • New or current test strip lots • New or current QC material L1 and L3 | <ul style="list-style-type: none"> • 50% Dextrose IV solution • Roche false bottom tubes | <ul style="list-style-type: none"> • Nova StatStrip Glucose Meter x2 • Roche Cobas c501 • EP Evaluator • Centrifuge • 37 ° Water or sand bath |

4.0 SPECIAL SAFETY PRECAUTIONS:

Exercise the normal precautions required for handling all laboratory reagents.
Disposal of all waste material should be in accordance with local guidelines.

5.0 PROCEDURE:

5.1. Validation of New Lot of Nova Glucose **Control Material 1 and 3** (day to day precision)

| Step | Action |
|------|--|
| 1 | Receive new lots of Control Levels 1 and 3 from Nova Biomedical. Prepare an excel spreadsheet similar to the one pictured in Appendix A. |
| 2 | Document the new lot numbers, expiry dates and expected ranges (on the front of the vials) on the spreadsheet. |
| 3 | Document the lot number and expiry date of the current test strip lot numbers used for testing. Prepare 2 Nova StatStrip Glucose Meters documenting the model numbers used. |
| 4 | Test the new levels 1 and 3 once a day on 2 separate meters for 10 days (total of 20 values for each level). Follow doc. #110-130-01 (StatStrip Blood Glucose POC Meter) for correct operation of the meter. QC values MUST fall within the ranges supplied. |
| 5 | Plot the data using the EP evaluator Simple Precision Module (see example Appendix B). File pdf in shared drive under "Point of Care Glucose Meters". |
| 6 | If all data fall within the acceptable range on the vials, the new QC lots are validated. If all data does not fall within range contact Clinical Biochemist for review and next steps. |

5.2. Validation of New Lot of **Glucose Test strips** (patient comparisons):

| Step | Action |
|------|--|
| 1 | Receive new lot of glucose test strips from Nova Biomedical. Prepare to do a 10 sample patient correlation (using the new strips) between the Nova StatStrip meter (in duplicate between 2 meters) and the Roche c501 (only run once). This usually takes a day to perform. |
| 2 | Obtain a 50% Dextrose IV solution from pharmacy. This will be used as a spiking solution. |
| 3 | Prepare an excel spreadsheet to enter all data similar to the one pictured in Appendix C. Document all lot numbers and expiration dates as listed. |
| 4 | On the morning of testing turn on the 37 °C heating bath. Run Daily QC on both meters using the new strips. |
| 5 | Fresh lithium heparin whole blood is required for the Nova glucose meter which is then immediately centrifuged to obtain plasma to analyze on the Roche c501. The sample values must span the meter's measuring range up to 30 mmol/L with 2 specific target values of 1.7 and 2.2 mmol/L. Fresh specimens can be gathered by removing 1 mL of whole blood from random samples arriving in the lab. Be sure to cap tightly when complete. |
| 6 | Aliquot the 1 mL whole blood specimen into a Roche false bottom tube and mix well. Quickly analyze and the whole blood on 1 meter to get a starting baseline value. |
| 7 | Specimens can be used as is OR are either spiked with the dextrose solution to increase the values as desired or glycolyzed in the heating bath to lower the values. Lower the values by placing the 1 mL whole blood sample into the 37 ° heating bath. The glucose value will decrease by 0.84 mmol/L for every hour it remains in the bath. Use this information to calculate the approximate amount of time you need to arrive at the target value. Refer to doc# 110-130-07 "Evaluation Protocol for Point of Care (POC) Blood Glucose Monitoring Instruments" section 7 and 8 for information on spiking and glycolyzing whole blood specimens. |
| 8 | Continue adding samples so that the entire range is covered. Whole blood samples should be analyzed as soon as possible and the specimens spun immediately for analysis on the c501. |
| 9 | Once all data is collected enter the values into the EP Evaluator software "Glucose POC Instrument Evaluation" statistical module. See Appendix D. |
| 10 | Forward completed data and EP evaluator pdf report to a Biochemist for sign-off. |
| 11 | Notify Nova Biomedical that strip lot is approved. |


Appendix A

**Nova Biomedical – New QC
Evaluation**

Site : _____
 Analyst: _____
 Start Date: _____
 Test strip lot: _____ exp _____
 Statstrip Meter #1 SN# _____
 Meter #2 SN# _____

| | Meter #1 | | Meter #2 |
|--|-------------|--------------|--------------|
| Nova Low Glucose Control (1) lot # xxxxxxxxxxxx exp xxxxx <u>Allowable Range</u> xx- xx mmol/L (printed on the vial) | Date | Value | Value |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| Nova High Glucose Control (3) lot # xxxxxxxxxxxx exp xxxxx <u>Allowable Range</u> xx- xx mmol/L (printed on the vial) | Date | Value | Value |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

APPENDIX B



| | Source | Lot Number | Expiration | Analytes |
|----------|-----------------|------------------|-------------|----------|
| Controls | Nova Biomedical | 0417342301 | 08 Dec 2020 | Glucose |
| | Nova Biomedical | 0417346303 | 12 Jun 2020 | Glucose |
| Reagents | Nova Biomedical | Strip 0317129309 | 30 Nov 2019 | Glucose |

Precision Data

| Index | Results | Index | Results | Index | Results | Index | Results |
|-------|---------|-------|---------|-------|---------|-------|---------|
| 1 | 2.9 | 6 | 3.1 | 11 | 3.4 | 16 | 3.3 |
| 2 | 3.4 | 7 | 3.2 | 12 | 2.8 | 17 | 3.3 |
| 3 | 3.3 | 8 | 3.2 | 13 | 3.4 | 18 | 3.1 |
| 4 | 3.3 | 9 | 3.3 | 14 | 3.0 | 19 | 3.3 |
| 5 | 3.1 | 10 | 3.5 | 15 | 2.7 | 20 | 3.2 |

O: outliers X: excluded from calculations

Appendix C

| Evaluation of Nova Test strips – Patient correlation vs Roche Cobas | | | |
|---|--|-----------|--|
| New test strip Lot # | | Exp. date | |
| Site: | | | |
| Analyst: | | | |
| Date of correlation: | | | |
| Meter Serials # | | | |
| Nova Control 1 lot # | | Exp. date | |
| Nova Control 3 lot # | | Exp. date | |
| 50% glucose lot # | | Exp. date | |
| Specimen type: | | | |

| | | Patient Correlation | | | | | |
|--------|------|----------------------|----------|----------|---|---------------------------|----------|
| Spec # | Date | Nova Stat Strip (WB) | | | | Roche Cobas c501 (plasma) | |
| | | Analysis time | Result 1 | Result 2 | X | Analysis time | Result 1 |
| 1 | | | | | | | |
| 2 | | | | | | | |
| 3 | | | | | | | |
| 4 | | | | | | | |
| 5 | | | | | | | |
| 6 | | | | | | | |
| 7 | | | | | | | |
| 8 | | | | | | | |
| 9 | | | | | | | |
| 10 | | | | | | | |
| 11 | | | | | | | |

mmol/L

APPENDIX D

