


Document History:

Title: StatStrip Glucometer Operator's Manual **Site(s):** Shared Health

Document #:	110-130-10	Version #:	01
Section:	Clinical Biochemistry	Subsection:	Point of Care

Approved by:	Hayden Malvern	Date:	06-June-2018
Signature:		Effective Date:	06-July-2018

Details of Recent Revision

1. New document.

Hayden Malvern 06-June 2018

Purpose	To provide direction on safe use of Blood Glucose Meters (Nova StatStrip), and associated infection control and quality assurance requirements for laboratory staff.
Introduction	Glucose is the major carbohydrate present in the peripheral blood. Oxidation of glucose is the major source of cellular energy in the body. Glucose derived from dietary sources is converted to glycogen for storage in the liver or to fatty acids for storage in adipose tissue. The concentration of glucose in blood is controlled within narrow limits by many hormones, the most important of which are produced by the pancreas. The most frequent cause of hyperglycemia is diabetes mellitus resulting from a deficiency in insulin secretion or action. A number of secondary factors also contribute to elevated blood glucose levels. These include pancreatitis, thyroid dysfunction, renal failure and liver disease.
Principle	<p>The glucose measurement is based on the following methodology:</p> $\begin{array}{l} \text{Glucose} + \text{Enzymes(oxidized form)} \rightarrow \text{Gluconic Acid} + \text{Enzymes(reduced form)} \\ \text{Enzymes(reduced form)} + \text{Ferricyanide} \rightarrow \text{Enzymes(oxidized form)} + \text{Ferrocyanide} \\ \text{Ferrocyanide} \xrightarrow[\text{Electrode}]{-e^-} \text{Ferricyanide} \end{array}$ <p>The current generated at the electrode is proportional to the glucose concentration of the sample.</p>
Definitions	<p>Blood Glucose Meter – A portable electronic device used to measure the capillary blood glucose level.</p> <p>Capillary Blood Glucose – The amount of glucose in blood is measured on a sample obtained by lancing the skin (usually the fingertip or alternate site, e.g. the palmar area of the thumb).</p>
Device Training	Staff using meter must receive training in proper use of the StatStrip Blood Glucose Meter, and associated quality assurance processes. Assign MediaLab TACA, 230-10-57.
Contraindications	Capillary Blood Glucose testing may not be clinically appropriate when peripheral blood flow is decreased. Shock, severe hypotension, hyperosmolar hyperglycemia and occurrence of severe dehydration are examples of clinical conditions that may adversely affect the measurement of glucose in the peripheral blood.
Use & Reporting	Laboratory staff use the StatStrip Blood Glucose Meter to check the patient's blood glucose concentration in order to determine if the result meets the criterion for safety for administering the glucose drink, required for performing the Oral Glucose Tolerance Test (OGTT).

Storage & Handling	Blood Glucose Meters, test strips and quality control solutions must be stored according to manufacturer's instructions.		
Materials		Vendor code	SAP
	StatStrip Glucose Meter	53400	205978
	Glucose test strip	42214	203423
	Alcohol swabs		
	Lancets		
	Control Solutions L1/ L3	41741, 41743	205980, 205979
	StatStrip Glucose Linearity Kit	42173	302261
Safety	Follow Routine Laboratory practice		
Sample Requirement	Whole Blood (Fresh, Sodium heparin, Lithium heparin): Arterial, Venous, Capillary		
Quality Control	<p>Level 1 and Level 3 of Nova StatStrip Glucose control solutions are run and evaluated once every 24 hours or before patient or EPT testing whichever is less frequent; and if a problem is suspected with the meter.</p> <p>EPT – Subscribe to CAP- Quality Cross Check -Whole Blood Glucose (WBGQ) program</p>		
Linearity Check	HSC & SBH Labs Only – Run the linearity check solutions every 6 months. Evaluate against acceptance criteria and document.		
Running Samples	Follow the procedure below to run and evaluate the quality control solutions		
	Step		
	1	From the Patient Test screen	
		If	Then
		QC	press the QC
		Patient sample	press the Accept
	2	The Enter Strip Lot screen is displayed. Enter the Strip Lot Number or scan the barcode. Press the Accept, if the lot number is correct.	
	3	If	Then
		QC	The Enter QC Lot screen is displayed. Enter the QC lot number or scan the barcode. Press the Accept, if the lot number is correct.
		Patient sample	From the Enter Patient ID screen, enter or scan the Patient ID. Once the Patient's ID/Accession Number has been entered, press the Accept.
	4	The Insert Strip screen is displayed. Insert a Test Strip	

	5	With the test strip correctly inserted, the Apply Sample screen is displayed.	
	6	If	Then
		QC	-Gently mix the StatStrip Glucose Control Solution. -Discard the first drop of control solution from the bottle to avoid contamination. -Place a drop of control solution from the bottle at the end of the test strip until the well of the test strip is full and the meter beeps -Recap the control solution.
		Patient sample	Following SOP 100-10-79 "Phlebotomy Collection Manual", collect a skin puncture sample . Wipe away the first drop of a skin puncture sample. Touch the end of the test strip to the blood drop until the well of the test strip is full and the meter beeps.
	<p>CAUTION: If the test strip does not fill completely, DO NOT touch the test strip to the blood a second time. Discard the strip and start over with a new strip.</p>		
	7	The Testing Sample screen is displayed. In six seconds, the test results are displayed in mmol/L.	
	8	Use the ejector button on the back of the meter to eject the strip directly into a biohazard container.	
	9	To accept the result, press Accept.	
	10	Evaluate and document results	
		If	Then
		QC results	Evaluate and document in QC log
		Patient results	Slot enter the result in Delphic LIS
	11	When the meter is not in use, place it into the Docking/Charging Station.	
Reporting		Slot enter results into the LIS	
Reference Range		3.3 – 6.0 mmol/L	
Reportable Range		0.6 - 33.3 mmol/L	

Clinical Limits	WRHA	Critical Low	Normal Low	Normal High	Critical High
	Adult	2.5	3.6	6.0	25.0
	Pediatric	2.5	3.3	5.6	25.0
	Neonate	1.7	2.8	4.4	18.0
	<p>Rational: These limits have been programmed into all Glucose meters by Biomedical Engineering for clinical use. Meters will flag when limits are breached.</p> <p>Action: Draw a random or FBS sample, accession, analyze on the I-stat or Automated lab analyzer and report as random or FBS as applicable.</p>				
Maintenance	As Required				
Replacing the Battery	Step				
	1	Press the Power button to enter the Sleep Mode.			
		If	Then		
		Replace within 20 seconds	Date/Time settings are not lost		
		>20 seconds	Refer to operators manual to set Date/Time		
	2	Pull back on the cover latch to release the cover. Take the battery cover off the back of the meter.			
	3	Push up on the battery latch. Remove the drained battery.			
	4	Replace with a fully charged battery.			
	5	Replace the battery cover.			
	6	Place the drained battery into the Charging Station.			
Cleaning the Meter		Apply 70% Isopropyl Alcohol to a soft cloth and wipe the meter surface. Once complete, immediately dry thoroughly.			

<p>Interfering Substances</p>	<p>The StatStrip Glucose Hospital Meter exhibits no interference from the following substances up to the following concentration levels:</p>																																																																					
	<table border="1"> <thead> <tr> <th data-bbox="389 394 738 430">Interfering Substances</th> <th colspan="2" data-bbox="747 394 1344 430">Concentration Level</th> </tr> </thead> <tbody> <tr> <td data-bbox="389 430 738 466">Acetaminophen</td> <td data-bbox="747 430 1096 466">10.0 mg/dL</td> <td data-bbox="1104 430 1344 466">0.66 mmol/L</td> </tr> <tr> <td data-bbox="389 466 738 501">Ascorbic Acid</td> <td data-bbox="747 466 1096 501">10.0 mg/dL</td> <td data-bbox="1104 466 1344 501">0.57 mmol/L</td> </tr> <tr> <td data-bbox="389 501 738 537">Bilirubin</td> <td data-bbox="747 501 1096 537">15.0 mg/dL</td> <td data-bbox="1104 501 1344 537">0.26 mmol/L</td> </tr> <tr> <td data-bbox="389 537 738 573">Cholesterol</td> <td data-bbox="747 537 1096 573">500.0 mg/dL</td> <td data-bbox="1104 537 1344 573">12.9 mmol/L</td> </tr> <tr> <td data-bbox="389 573 738 609">Creatinine</td> <td data-bbox="747 573 1096 609">6.0 mg/dL</td> <td data-bbox="1104 573 1344 609">0.53 mmol/L</td> </tr> <tr> <td data-bbox="389 609 738 644">Dopamine</td> <td data-bbox="747 609 1096 644">10.0 mg/dL</td> <td data-bbox="1104 609 1344 644">0.53 mmol/L</td> </tr> <tr> <td data-bbox="389 644 738 680">Ephedrine</td> <td data-bbox="747 644 1096 680">0.9 mg/dL</td> <td data-bbox="1104 644 1344 680">0.055 mmol/L</td> </tr> <tr> <td data-bbox="389 680 738 716">D(+) Galactose</td> <td data-bbox="747 680 1096 716">350.0 mg/dL</td> <td data-bbox="1104 680 1344 716">19.4 mmol/L</td> </tr> <tr> <td data-bbox="389 716 738 751">Hematocrit (RBC)</td> <td data-bbox="747 716 1096 751">20% - 65%</td> <td data-bbox="1104 716 1344 751"></td> </tr> <tr> <td data-bbox="389 751 738 787">Ibuprofen</td> <td data-bbox="747 751 1096 787">48.0 mg/dL</td> <td data-bbox="1104 751 1344 787">2.33 mmol/L</td> </tr> <tr> <td data-bbox="389 787 738 823">L-Dopa</td> <td data-bbox="747 787 1096 823">100.0 mg/dL</td> <td data-bbox="1104 787 1344 823">5.07 mmol/L</td> </tr> <tr> <td data-bbox="389 823 738 900">D(+) Maltose Monohydrate</td> <td data-bbox="747 823 1096 900">240.0 mg/dL</td> <td data-bbox="1104 823 1344 900">6.66 mmol/L</td> </tr> <tr> <td data-bbox="389 900 738 936">D(+) Maltotetraose</td> <td data-bbox="747 900 1096 936">240.0 mg/dL</td> <td data-bbox="1104 900 1344 936">3.6 mmol/L</td> </tr> <tr> <td data-bbox="389 936 738 972">D(+) Maltotriose</td> <td data-bbox="747 936 1096 972">240.0 mg/dL</td> <td data-bbox="1104 936 1344 972">4.76 mmol/L</td> </tr> <tr> <td data-bbox="389 972 738 1008">Methyl-Dopa</td> <td data-bbox="747 972 1096 1008">1.0 mg/dL</td> <td data-bbox="1104 972 1344 1008">0.042 mmol/L</td> </tr> <tr> <td data-bbox="389 1008 738 1043">Oxygen</td> <td colspan="2" data-bbox="747 1008 1344 1043">All Concentrations</td> </tr> <tr> <td data-bbox="389 1043 738 1079">Salicylate</td> <td data-bbox="747 1043 1096 1079">30.0 mg/dL</td> <td data-bbox="1104 1043 1344 1079">1.87 mmol/L</td> </tr> <tr> <td data-bbox="389 1079 738 1115">Tetracycline</td> <td data-bbox="747 1079 1096 1115">30.0 mg/dL</td> <td data-bbox="1104 1079 1344 1115">0.62 mmol/L</td> </tr> <tr> <td data-bbox="389 1115 738 1150">Tolazamide</td> <td data-bbox="747 1115 1096 1150">15.0 mg/dL</td> <td data-bbox="1104 1115 1344 1150">0.48 mmol/L</td> </tr> <tr> <td data-bbox="389 1150 738 1186">Tolbutamide</td> <td data-bbox="747 1150 1096 1186">45.0 mg/dL</td> <td data-bbox="1104 1150 1344 1186">1.67 mmol/L</td> </tr> <tr> <td data-bbox="389 1186 738 1222">Triglycerides</td> <td data-bbox="747 1186 1096 1222">750.0 mg/dL</td> <td data-bbox="1104 1186 1344 1222">8.78 mmol/L</td> </tr> <tr> <td data-bbox="389 1222 738 1257">Uric Acid</td> <td data-bbox="747 1222 1096 1257">20.0 mg/dL</td> <td data-bbox="1104 1222 1344 1257">1.05 mmol/L</td> </tr> </tbody> </table>		Interfering Substances	Concentration Level		Acetaminophen	10.0 mg/dL	0.66 mmol/L	Ascorbic Acid	10.0 mg/dL	0.57 mmol/L	Bilirubin	15.0 mg/dL	0.26 mmol/L	Cholesterol	500.0 mg/dL	12.9 mmol/L	Creatinine	6.0 mg/dL	0.53 mmol/L	Dopamine	10.0 mg/dL	0.53 mmol/L	Ephedrine	0.9 mg/dL	0.055 mmol/L	D(+) Galactose	350.0 mg/dL	19.4 mmol/L	Hematocrit (RBC)	20% - 65%		Ibuprofen	48.0 mg/dL	2.33 mmol/L	L-Dopa	100.0 mg/dL	5.07 mmol/L	D(+) Maltose Monohydrate	240.0 mg/dL	6.66 mmol/L	D(+) Maltotetraose	240.0 mg/dL	3.6 mmol/L	D(+) Maltotriose	240.0 mg/dL	4.76 mmol/L	Methyl-Dopa	1.0 mg/dL	0.042 mmol/L	Oxygen	All Concentrations		Salicylate	30.0 mg/dL	1.87 mmol/L	Tetracycline	30.0 mg/dL	0.62 mmol/L	Tolazamide	15.0 mg/dL	0.48 mmol/L	Tolbutamide	45.0 mg/dL	1.67 mmol/L	Triglycerides	750.0 mg/dL	8.78 mmol/L	Uric Acid	20.0 mg/dL
Interfering Substances	Concentration Level																																																																					
Acetaminophen	10.0 mg/dL	0.66 mmol/L																																																																				
Ascorbic Acid	10.0 mg/dL	0.57 mmol/L																																																																				
Bilirubin	15.0 mg/dL	0.26 mmol/L																																																																				
Cholesterol	500.0 mg/dL	12.9 mmol/L																																																																				
Creatinine	6.0 mg/dL	0.53 mmol/L																																																																				
Dopamine	10.0 mg/dL	0.53 mmol/L																																																																				
Ephedrine	0.9 mg/dL	0.055 mmol/L																																																																				
D(+) Galactose	350.0 mg/dL	19.4 mmol/L																																																																				
Hematocrit (RBC)	20% - 65%																																																																					
Ibuprofen	48.0 mg/dL	2.33 mmol/L																																																																				
L-Dopa	100.0 mg/dL	5.07 mmol/L																																																																				
D(+) Maltose Monohydrate	240.0 mg/dL	6.66 mmol/L																																																																				
D(+) Maltotetraose	240.0 mg/dL	3.6 mmol/L																																																																				
D(+) Maltotriose	240.0 mg/dL	4.76 mmol/L																																																																				
Methyl-Dopa	1.0 mg/dL	0.042 mmol/L																																																																				
Oxygen	All Concentrations																																																																					
Salicylate	30.0 mg/dL	1.87 mmol/L																																																																				
Tetracycline	30.0 mg/dL	0.62 mmol/L																																																																				
Tolazamide	15.0 mg/dL	0.48 mmol/L																																																																				
Tolbutamide	45.0 mg/dL	1.67 mmol/L																																																																				
Triglycerides	750.0 mg/dL	8.78 mmol/L																																																																				
Uric Acid	20.0 mg/dL	1.05 mmol/L																																																																				
<p>New or Replacement Meters</p>	<p>Prior implementing for patient testing.</p>																																																																					
	<table border="1"> <thead> <tr> <th data-bbox="389 1407 868 1470">HSC & SBH</th> </tr> </thead> <tbody> <tr> <td data-bbox="389 1470 868 1911"> <p>Run all 5 levels of the Linearity solutions and the two level of QC. Evaluate against acceptance criteria and document for accreditation purposes</p> </td> </tr> </tbody> </table>	HSC & SBH	<p>Run all 5 levels of the Linearity solutions and the two level of QC. Evaluate against acceptance criteria and document for accreditation purposes</p>	<table border="1"> <thead> <tr> <th data-bbox="885 1407 1344 1470">All Other Sites</th> </tr> </thead> <tbody> <tr> <td data-bbox="885 1470 1344 1911"> <p>Request the validation documentation from your Biomedical Engineering when you pick up the meter.</p> <p>Run the two levels of QC. Evaluate against acceptance criteria and document with the validation documentation from Biomedical Engineering.</p> </td> </tr> </tbody> </table>	All Other Sites	<p>Request the validation documentation from your Biomedical Engineering when you pick up the meter.</p> <p>Run the two levels of QC. Evaluate against acceptance criteria and document with the validation documentation from Biomedical Engineering.</p>																																																																
HSC & SBH																																																																						
<p>Run all 5 levels of the Linearity solutions and the two level of QC. Evaluate against acceptance criteria and document for accreditation purposes</p>																																																																						
All Other Sites																																																																						
<p>Request the validation documentation from your Biomedical Engineering when you pick up the meter.</p> <p>Run the two levels of QC. Evaluate against acceptance criteria and document with the validation documentation from Biomedical Engineering.</p>																																																																						


Troubleshooting	Refer to the Nova Instructions for Use Manual Call clinical Biomedical department for your site or region Contact Nova Biomedical Technical Support toll free at 1-800-545-NOVA.
Associated Documents	MediaLab TACA, 230-10-57 QC for StatStrip Blood Glucose POC Meter - F110-10-17
Reference	StatStrip Glu IFU 1.86 PN53083C EN DSM Lab Information Manual SOP 110-130-01 v2 WRHA Program clinical ranges (HSC Biomedical Engineering)

StatStrip®

Glucose Monitoring System Quick Operating Guide




1 From Home screen, press Login.




2 Enter or scan Operator ID and press Accept.




3 From Patient Test screen, press Accept.




4 Check Strip lot no. and press Accept.



5 Enter or scan Patient ID and press Accept.




6 Insert Test Strip into Meter.



7 Wash patient's hand thoroughly and massage finger to stimulate blood flow.




8 Use safety lancet to puncture finger / squeeze finger to form blood drop. Wipe away the first blood drop.



9 Touch strip to blood drop. Result will appear within 8 seconds.

10 **Warning!**

The test strip must fill completely upon touching the blood droplet. If the test strip does not fill completely, do not touch the blood droplet a second time. Discard the test strip and repeat the test with a new test strip.



11 To accept result, press Accept. To reject result, press Reject.



12 To review other results, press Review from Patient Test screen.




Nova Biomedical Canada Ltd., 17 - 2900 Argentia Road,
 Mississauga, ON L5N 7X9 Canada
 Tel: 1-800-263-5999 • www.novabiomedical.com

StatStrip®

Glucose Monitoring System Quick QC Guide




1 From Home screen, press Login.




2 Enter or scan Operator ID and press Accept.




3 From Patient Test screen, press QC.



4 Check Strip lot no. and press Accept.



5 Enter QC lot no. and press Accept.




6 Insert Test Strip into Meter.




7 Touch drop from QC bottle to strip. Result will appear within 6 seconds.

8 **Warning!**

The test strip must fill completely upon touching the QC droplet. Do not add a second QC drop to the test strip. Discard the test strip and repeat the test with a new test strip.



9 Result will appear within 6 seconds.




10 To accept result, press Accept.




Nova Biomedical Canada Ltd., 17 - 2900 Argentinia Road,
 Mississauga, ON L5N 7X9 Canada
 Tel: 1-800-263-5999 • www.novabiomedical.com

StatStrip®


Glucose Monitoring System Troubleshooting Guide




1 Change battery or place meter into charging/docking station.



2 The test has been cancelled, repeat the test with a new test strip. Leave the test strip in place until the result is displayed on the screen.




3 Meter will only work in a temperature range of 59-104°F (15-40°C).



4 Insert new strip and repeat test. If the error code persists, perform the test using an alternate strip vial or alternate method.




5 Occurs after test strip insertion or during analysis. Insert another strip and repeat the test. If the error code persists, perform the test using an alternate strip vial or alternate method.



6 Either insufficient sample or the sample was applied incorrectly. Repeat test with a new strip. If the error code persists, perform the test using an alternate method.



7 Server refuses to allow dialog with meter or connection to server was broken. Check network settings, network status, or contact your administrator for assistance.



8 Meter was removed before data transfer was complete. Re-dock the meter.



Nova Biomedical Canada Ltd., 17 - 2900 Argentia Road,
 Mississauga, ON L5N 7X9 Canada
 Tel: 1-800-263-5999 • www.novabiomedical.com

APPENDIX A:

Stat Strip Glucose Monitoring System
Quick User Guide

<u>Maintenance</u>	<u>Quality control (QC)</u>	<u>Patient Testing</u>
<ol style="list-style-type: none"> 1. Clean Glucose meter between Patient use with OxyvirTB wipes. 2. Dock meter in docking station when not in use. 3. Test Strips expire at 6 months once vial open. Date vial once opened. 4. QC bottles expire at 3 months once vial opened. Date vial once opened. 5. When testing, keep meter flat. Remove test stripe before accepting results. <i>Tilting meter up may result in fluid entering and damaging the meter.</i> 	<ol style="list-style-type: none"> 1. QC to be completed every 24 Hrs. Visual message warning when QC due. 2. From welcome screen, press login. 3. Enter operator Id (1) & press Accept. 4. From Pt test screen, press QC. 5. Scan test stripe lot no. 6. Scan QC lot no. 7. Insert test strip into meter, touch drop from QC bottle to tip of strip. QC result within range, press accept. 8. Complete Qc with both Qc bottles (Low and High) 	<ol style="list-style-type: none"> 1. From welcome screen, press login. 2. Enter operator Id (1) & press accept. 3. From patient test screen, press Accept. 4. Scan test stripe lot no. 5. Enter Patient Id (1) , press Accept. 6. Insert test strip into meter, touch tip of strip to pt blood drop. Do not add second Blood drop. Keep meter flat or tilted downward. 7. Wait 6 sec. Choose Accept or reject test results. 8. Meter range allow low-0.5 , high-33mmol/l.

M. Ash CRN ICMS Nov5/10

APPENDIX B:**Running a Linearity Test**

1. From the Patient Test screen, press Menu.
2. From the Menu screen, press Linearity.
3. The Enter Strip Lot screen displays. Enter the Strip Lot Number or scan the barcode.
4. Press Accept if the lot number is correct.
5. The Enter Linearity Lot screen displays. Enter the Linearity lot number, select from the Linearity Lot List screen (press the List button), or scan the barcode.
6. Press Accept if the lot number is correct.
7. The Insert Strip screen displays. Insert a Test Strip as shown on the screen.
8. With the test strip correctly inserted, the Apply Sample screen displays.
9. Gently mix the StatStrip Linearity Solution before each use.
10. Discard the first drop of linearity solution from the bottle to avoid contamination.
11. Place a drop of linearity solution from the bottle at the end of the test strip until the solution is drawn into the well of the test strip. When enough sample has been drawn into the strip, an audible beep is sounded by the meter.
12. Recap the linearity solution. The Testing Sample screen displays. The screen shows a clock with seconds remaining below the clock.
13. When the meter completes the test, the Linearity Result screen displays with the results.
14. Document the result on the Linearity Check & Meter validation Data Log.
15. Remove the strip manually or use the ejector button on the back of the meter to eject the strip directly into a biohazard container.