

Year

Month

## **Clinitek Status Maintenance Log**

Document # F110-50-19
Version # 02

Approved by:											
Haydullalven											
H Malvern											

Site

Effective Date: 21 May 2020

Serial #

Source Document: 110-50-19

## Laboratory Maintenance Record Clinitek Status Urinalysis Analyzer

Daily	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Clean Feed table																															
Clean Surface																															
Perform QC																															
Tech's initials																															
																1		1		•											
Weekly								Week 1 Date					Week 2 Date					Week 3 Date					Week 4 Date					Week 5 Date			
Clean White Calibra																															
Tech's Initials																															
												Date	)									Tech's Initials / Comments									
Unscheduled Maintenance / As required																															
Changing the Paper																															
Disinfect the Instrument																															
											,										,										
Reviewed by												_ Da	ate_																		