

ERYTHROCYTE SEDIMENTATION RATE (ESR) ANALYZER SPECIMENS

Date Issued: March 15, 2019

Westman Laboratory will go live with the new analyzer on Saturday, March 16, 2019. Health Sciences Centre site is still completing the validation.

The sample required can be a 2, 3, or 4 mL EDTA tube, but must have a minimum volume of 2 mL.

The stability of the sample has not changed
ie. 36 hours for refrigerated samples at 2 – 8°C or 24 hours for room temperature samples at 18 – 25°C.

Please contact Dr. Ping Sun at 204-787-4682 or psun@sharedhealthmb.ca or
Cathy McMahon, 204-787-4971 or cmcmahon@sharedhealthmb.ca with any questions or concerns.