

## HEPARIN INDUCED THROMBOCYTOPENIA (HIT) TESTING

### CHANGE IN SPECIMEN REQUIREMENT

**Date Effective: Immediately**

**Date Issued: July 4, 2018**

#### **Change in or New Test Procedure:**

The McMaster Platelet Immunology Laboratory has made changes in their HIT testing protocol and sample requirements. Therefore, additional collection of blood/tubes will be needed, should the result of our local test be positive or indeterminate.

**ALL** the following tubes **must** be obtained for comprehensive testing:

- **Serum: 1 red top 10 mL tube (at least 2/3 full)**
- **Plasma: 3 light blue 1.8 mL Sodium Citrate tubes**

These requirements will be re-assessed once the new local HIT testing method is implemented.

Shared Health testing site: SBH Hematology Laboratory

Please use requisition below.

[H.I.T. \(Heparin-Induced Thrombocytopenia\) Screening Requisition F140-140-26A](#)

#### **More information:**

See the complete details of the Heparin Induced Thrombocytopenia (HIT) in the LAB Information Manual (LIM) below:

[HIT](#)

All Shared Health Clinical Practice Changes can be found at:

[Clinical Practice Changes](#)

#### **Contact Information:**

Carmen Morales, [cmorales@sharedhealthmb.ca](mailto:cmorales@sharedhealthmb.ca) Phone: 204-787-4682

-or-

Cathy McMahan, [cmcmahan@sharedhealthmb.ca](mailto:cmcmahan@sharedhealthmb.ca) Phone: 204-787-4971