



Document #: CO-4	 Standard Operating Procedure (SOP)	Status: v 0.0 Revised: June 01, 2020 Owner: Shared Health
Reprocessing Single Use N95 Respirators due to COVID-19 using N95 Respirator Content Information from Vendors		Page 1 of 4

Purpose	To provide information from vendors on the cellulose, latex / natural rubber content for N95 respirators.
----------------	---


*Highlighted models are currently available in Manitoba as of date of document.

3M Based on 3M Technical Bulletin May, 2020 Revision 5 Cellulose Certification - Filtering Facepiece Respirators https://multimedia.3m.com/mws/media/18246130/cellulose-certification-filtering-facepiece-respirators.pdf Note - There is no latex in any 3M disposable respirator.	
The products listed immediately below contain cellulose.	Reprocessing Technology
<ul style="list-style-type: none"> • 3M™ VFlex™ Healthcare Particulate Respirator and Surgical Mask 1804, N95 • 3M™ VFlex™ Healthcare Particulate Respirator and Surgical Mask 1804S, N95, Small • 3M™ VFlex™ Health Care Particulate Respirator and Surgical Mask 1805, N95 • 3M™ VFlex™ Health Care Particulate Respirator and Surgical Mask 1805S, N95, Small • 3M™ Particulate Respirator 8000, N95 • 3M™ Particulate Respirator 8200/07023(AAD), N95 • 3M™ Particulate Welding Respirator 8212, N95 with Faceseal • 3M™ Particulate Respirator 8214, N95, with Faceseal and Nuisance Level Organic Vapor Relief • 3M™ Particulate Respirator 8233, N100 • 3M™ Particulate Respirator 8293, P100 • 3M™ Particulate Welding Respirator 8512, N95 • 3M™ Particulate Respirator 8514, N95, with Nuisance Level Organic Vapor Relief • 3M™ Particulate Respirator 8822, FFP2 • 3M™ Particulate Respirator 8822KV, P2 • 3M™ Particulate Respirator 8822AUS, P2 • 3M™ Particulate Respirator 9001, KN90 • 3M™ Particulate Respirator 9002, KN90 • 3M™ Particulate Respirator 9010, N95 • 3M™ Particulate Respirator 9010CN, N95 • 3M™ Particulate Respirator 9010CN, KN90 • 3M™ Particulate Respirator 9010MX, N95 • 3M™ Particulate Respirator 9010V, N95 • 3M™ VFlex™ Particulate Respirator 9105, N95 • 3M™ VFlex™ Particulate Respirator 9105S, N95 • 3M™ Particulate Respirator & Surgical Mask 9132, N95 • 3M™ Particulate Respirator 9501, KN95 • 3M™ Particulate Respirator 9502, KN95 	None

Document #: CO-4	 Standard Operating Procedure (SOP)	Status: v 0.0 Revised: June 01, 2020 Owner: Shared Health
---------------------------------------	---	---

Reprocessing Single Use N95 Respirators due to COVID-19 using N95 Respirator Content Information from Vendors	Page 2 of 4
--	-------------

The products listed immediately below have not been formulated to contain, and do not knowingly contain cellulose.	Reprocessing Technology
<ul style="list-style-type: none"> • 3M™ Health Care Particulate Respirator and Surgical Mask 1860, N95 • 3M™ Health Care Particulate Respirator and Surgical Mask 1860S, N95 • 3M™ Aura™ Health Care Particulate Respirator and Surgical Mask 1870, N95 • 3M™ Aura™ Health Care Particulate Respirator and Surgical Mask 1870+, N95 • 3M™ Particulate Respirator 8110S, N95 • 3M™ Particulate Respirator 8205J, DS2 • 3M™ Particulate Respirator 8210, N95 • 3M™ Particulate Respirator 8210CN, N95, KN95 • 3M™ Particulate Respirator 8210Plus, N95 • 3M™ Particulate Respirator 8210PlusCN, N95, KN95 • 3M™ Particulate Respirator 8210VCN, N95, KN95 • 3M™ Particulate Respirator 8210V, N95 • 3M™ Particulate Respirator 8211, N95 • 3M™ Particulate Respirator 8240, R95 • 3M™ Particulate Respirator 8246, R95, with Nuisance Level Acid Gas Relief • 3M™ Particulate Respirator 8247, R95, with Nuisance Level Organic Vapor Relief • 3M™ Particulate Respirator 8271, P95 • 3M™ Particulate Respirator 8510, N95 • 3M™ Particulate Respirator 8511, N95 • 3M™ Particulate Respirator 8511CN, N95 • 3M™ Particulate Respirator 8511CN, KN95 • 3M™ Particulate Welding Respirator 8515/07189(AAD), N95 • 3M™ Particulate Respirator 8516, N95, with Nuisance Level Acid Gas Relief • 3M™ Particulate Respirator 8576, P95, with Nuisance Level Acid Gas Relief • 3M™ Particulate Respirator 8577, P95, with Nuisance Level Organic Vapor Relief • 3M™ Particulate Respirator 8612F, N95 • 3M™ Particulate Respirator 8670F, N95 • 3M™ Disposable Cup Shape Healthcare Respirator 8801H, PFF-2(S) • 3M™ Particulate Respirator 8822KV, P2 • 3M™ Particulate Respirator 8822KV, P2 • 3M™ Aura™ Particulate Respirator 9205+, N95 • 3M™ Aura™ Particulate Respirator 9210/37021(AAD), N95 • 3M™ Aura™ Particulate Respirator 9210+/37192(AAD), N95 • 3M™ Aura™ Particulate Respirator 9211/37022(AAD), N95 • 3M™ Aura™ Particulate Respirator 9211+/37193(AAD), N95 • 3M™ Aura™ Particulate Respirator 9320+, FFP2, P2 • 3M™ Disposable Horizontal Flat Fold Respirator 9320+BR, PFF-2(S) • 3M™ Aura™ Particulate Respirator 9322+, FFP2 • 3M™ Particulate Respirator 9501+, KN95 • 3M™ Particulate Respirator 9501V+, KN95 • 3M™ Particulate Respirator 9502+, KN95 • 3M™ Particulate Respirator 9502+, N95 	Sterrad Stryker Steris CleanWorks


Document #: CO-4	 Standard Operating Procedure (SOP)	Status: v 0.0 Revised: June 01, 2020 Owner: Shared Health
---------------------------------------	---	---

Reprocessing Single Use N95 Respirators due to COVID-19 using N95 Respirator Content Information from Vendors	Page 3 of 4
--	-------------

<ul style="list-style-type: none"> • 3M™ Particulate Respirator 9502V+, KN95 • 3M™ Particulate Respirator 9505+, KN95 • 3M™ Particulate Respirator 9541, KN95 • 3M™ Particulate Respirator 9541V, KN95 • 3M™ Particulate Respirator 9542, KN95 • 3M™ Particulate Respirator 9542V, KN95 • 3M™ Particulate Respirator 9552, N95 • 3M™ Particulate Respirator 9552, KN95 • 3M™ Particulate Respirator 9552V, KN95 • 3M™ Disposable Vertical Flat Fold Respirator 9820 BR, PFF-2(S) • 3M™ Disposable Vertical Flat Fold Respirator 9820+BR, PFF-2(S) • 3M™ Disposable Vertical Flat Fold Healthcare Respirator 9920H, PFF-2(S) 	
---	--

Kimberly Clark	
The products listed immediately below are not tested for the presence of cellulose, cellulose is not intentionally added during the manufacturing or packaging process.	Reprocessing Technology
Kimberly-Clark Fluidshield N95 Particulate Filter Respirator and Surgical masks <ul style="list-style-type: none"> • KC46727 • KC46767 • KC46827 • KC46867 • KC62126 <p>Note - There is no latex in the above respirators</p>	Sterrad Stryker Steris Clean Works
Replacement elastics used for the Kimberly-Clark Fluidshield N95 Particulate Filter Respirator and Surgical masks contain latex/natural rubber.	Reprocessing Technology
<ul style="list-style-type: none"> • KC46827 MOD 	Sterrad Steris

Moldex	
Filtering facepiece respirators by Moldex do not contain cellulose.	Reprocessing Technology
<ul style="list-style-type: none"> • 1511 <p>Note - There is no latex in the above respirators</p>	Sterrad Stryker Steris Clean Works

Document #: CO-4	 Standard Operating Procedure (SOP)	Status: v 0.0 Revised: June 01, 2020 Owner: Shared Health
Reprocessing Single Use N95 Respirators due to COVID-19 using N95 Respirator Content Information from Vendors		Page 4 of 4

Pleats Plus by 3M	Reprocessing Technology
3M has been unable to verify cellulose content after purchase of Pleats Plus business. Clean Works has tested and indicated that these respirators can be reprocessed with the Clēan Flow Health Care Mini.	
<ul style="list-style-type: none"> • 50450 : 1054 • 50451: 1054S 	Clean Works

North by Honeywell	Reprocessing Technology
Vendor will not disclose content of masks	
<ul style="list-style-type: none"> • NON7130N95S • NON7130N95 	None

Dynamics	Reprocessing Technology
Vendor will not disclose content of masks	
<ul style="list-style-type: none"> • RPD714N95 	None