





SWIM – Access to Care

TWAIN Driver Install

Quick Reference Guide

Download and install the TWAIN driver needed for attaching and scanning documents within Novari ATC. This procedure only applies to workstations that are **not** managed by Shared Health. The driver will be installed automatically on managed Shared Health workstations.


| Install TWAIN Driver | |
|--|---|
| <p>1. Validate scanning support</p> | <ul style="list-style-type: none"> Log in to Novari ATC The application will launch and will automatically check for the TWAIN Driver on the workstation  |
| <p>2. Download the file</p> | <p>If the application determines that the driver is not installed, a pop-up will appear to download the driver</p> <ul style="list-style-type: none"> Click the Download button  |
| <p>3. Monitor download progress</p> | <p>The download progress will be displayed in the browser; if using Chrome it will display in the top right corner</p>  |
| <p>4. Start driver installation</p> | <ul style="list-style-type: none"> Once the download is complete, click the downloaded file to start the installation  |

Quick Reference Guide

Install TWAIN Driver


5. Complete driver installation

- Click the **Next** button on the **Setup Wizard** pop-up and follow the onscreen prompts
- Click **Close** when the installation is complete



6. Re-scan for TWAIN driver

- After the install window closes, refresh your browser or select the **“here”** link to re-initiate the TWAIN driver check



If the installation was successful the pop-up message will no longer appear when starting the application.

If the pop-up re-appears, repeat steps 2 – 6 or contact the Shared Health Service Desk



If you change workstations, you may be prompted to install the TWAIN driver on the new workstation.