


# ORDER OF DRAW CHART

## INVERSIONS:

Tubes are pictured in order of draw.



Red/Clear    Blue/Clear    Red    Gold    Red/Grey    Dark Green    Light Green    Lavender    Dark Blue    Yellow    Grey



## KNOW YOUR FILL LEVEL

YOUR FILL LEVEL HAS +/- 10% TOLERANCE  
QUALITY RESULTS = FULL TUBES  
ANYTHING ABOVE OR BELOW THE TOLERANCE LEVEL CAN IMPACT PATIENT TEST RESULTS

F100-10-20 V4.1

Effective date: 20 Oct 2022

Source document: 100-10-79