Physician Alert







May 20, 2015

Complement Activity, Total (CH50)

Effective date: May 25, 2015

The following changes will apply for <u>Complement Activity</u>, <u>Total (CH50)</u> measurement:

- Methodology from current EIA to <u>turbidimetry</u>
- Reference range 41.68 95.06 U/ml (adult). Range not established for pediatric age groups
- Availability: Monday Friday 8:00 a.m. 4:00 p.m.
- Average turn-around-time: 1 week
- Sample requirements: Serum 1mL minimum. <u>Must be separated and</u>
 frozen within 1 hour of collection. Please see instructions in the DSM LIM
- Performing site: <u>HSC Immunology</u>

Any questions or concerns may be directed to:

Dr. Carmen Morales, Medical Director, Hematology and Immunology Disciplines, 204-787-4682, cmarrier description descript

Sheila Ozamoto, Technical Director, Immunology, 204-787-3407, sozamoto@dsmanitoba.ca